


## Practical activities - Laboratory work: Osmosis in an onion – Allium cepa cells

|   |   |                    |
|---|---|--------------------|
| Student name:   | Class:  | Year:<br>2018/2019 |
| <b><u>Name of the project: CONSERVATION OF FOOD</u></b> |   |                    |
| <b>Theme:</b>   | <b>Osmosis in an onion – Allium cepa cells</b>  |                    |
| <b>Task:</b>  | Observe the influence of osmosis relation to an onion cells   |                    |
| <b>Materials and tools:</b>                             | onion skin – red onion epidermis, microscope, microscopic tools, slide and cover slip, paper towels, dropper, knife, 2 beakers, sodium chloride, (salt), water  |                    |
| <b>Working process:</b>                                 | <p><u>Observe the influence of osmosis relation to an onion cells</u></p> <ol style="list-style-type: none"> <li>1. Pell off a very, very thin piece of red onion using preparation needle. Make a wet mount of the red onion epidermis. Examine under microscope.</li> <li>2. After that take a dropper and add several drops of salt solution NaCl to one side of your cover slip, then take a small piece of paper towel and place it along the opposite edge of the cover slip. The paper should draw out the water and draw in the salt solution.</li> <li>3. Observed the effects of the saline (salt) solution on the onion cells about 5minutes.</li> <li>4. Make a drawing of the observed results.</li> </ol> |                    |

**Development  
off method::**


Observe the influence of osmosis relation to an onion cells

*a) one cell of an onion skin*



*magnification 10 x 20*

*b) one cell of the saline (salt) solution NaCl*



*magnification 10 x 20*

**Conclusion:**

- *In a concentrated salt solution NaCl the cell loses water and goes flaccid because the vacuole becomes flaccid and the cytoplasm stops pushing against the cell wall. This state is call **plasmolysis**.*

|  |  |
|--|--|
|  |  |
|  |  |